



made in the USA
with US & imported parts



15-4288 & 15-4288QS PHYSICAL THERAPY TABLE


User Guide



Shown in Pewter fabric




Chair
Position



Backrest

0° to 85°



Head Section
+40° to -90°

For complete product details please see User Manual

PRODUCT USE DESCRIPTION / IMPORTANT SAFETY INSTRUCTIONS

PRODUCT USE DESCRIPTION

The CanDo PT Tables are used to support the patient during examinations and procedures. They are intended to be operated by a healthcare professional in a medical environment. No special training is required but a review of the following Safety Instructions is important for the safety of the operator and patient. The healthcare professional should read and understand this entire manual before use with a patient. There are no known contraindications to the use of this equipment.

SYMBOL IDENTIFICATION



This symbol, when used in this manual and on product labels, represents a caution warning. Be sure to read and comply with all precautions and warnings.



This symbol, when used in this manual and on product labels, warns against an electrical shock hazard. Be sure to observe and comply with all warnings.



This symbol, when used in this manual or on product labels, indicates a Protective Earth (Ground) Terminal.



This symbol, when used in this manual and on product labels, indicates the name and address of the manufacturer.



This symbol, when used in this manual or on product labels, indicates the country of manufacture along with date of manufacture of the device next to it.



This symbol, when used in this manual and on product labels, indicates that the table and components are a Type B Applied Part pursuant to IEC 60601-1.



This symbol, when used in documentation on this product, indicates alternating current (AC).



This symbol, when used in documentation on this product, indicates direct current (DC).



This symbol is used to indicate that the operator should consult the user manual.

IMPORTANT SAFETY INSTRUCTIONS



CAUTION

READ AND SAVE THESE INSTRUCTIONS

The use of accessories, transducers, and cables other than those specified by the manufacturer, may result in increased emissions or decreased immunity of the table.

The Table should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the table should be observed to verify normal operation in the configuration in which it will be used.

The Table is designed to be a stand-alone table. This table must not be modified or incorporated into any other equipment.

As with any moving mechanism there are potential pinch points around and underneath the table. It is the responsibility of the operator of this equipment to insure that bystanders are not in the area below or around this equipment during operation.

IMPORTANT SAFETY INSTRUCTIONS

Proper operation of this equipment is very important for the safety of the operator, patient, and any other individuals in the area of this equipment. Directions for use of this equipment are described in this manual. The operator should read these sections carefully.

Weight Limit: (patient and accessories) 550 lbs. (250 kg) Do not exceed.

The Back Rest and Leg Rest sections are not designed to support the entire weight of the patient. Do not sit on the Back Rest or Leg rest sections.

Be certain that the table is completely lowered without any tilt being present prior to discharging an ambulatory patient. The patient may lose balance and fall.

This table is not suitable for use in the presence of flammable anesthetic mixture with air or oxygen or nitrous oxide.

Do not lift table by the table top. This can damage the table.

When lowering the table, make sure there is nothing underneath the table top that can impede motion (like stools, cabinets, accessory parts, cleaners, etc)

Use this furnishing only for its intended use as described in these instructions. Do not use attachments not recommended by the manufacturer.

Close supervision is necessary when this furnishing is used near children or disabled persons.



WARNING

To reduce the risk of burns, fire, electric shock or injury to persons:

1. Unplug this furnishing from the electrical outlet before cleaning.
2. Unplug from outlet before adding or removing parts.
3. Never operate this furnishing if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Contact CanDo Customer Service before use.
4. Keep the cord away from heated surfaces.
5. Never drop or insert any object into any opening.
6. Do not use outdoors.
7. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.



DANGER

Risk of electric shock - Connect this furnishing to a properly grounded outlet only. See Grounding Instructions in this manual.

Electrical Shock Hazard - The power supply/control module is located under the table seat. No user serviceable parts are inside the control box. Refer servicing to qualified personnel. Unplug wall plug prior to contact with any cables connected to the power supply.

INSTALLATION

The Tables come completely assembled and ready to use. Plug the cord into a functioning outlet that is rated for the table. (See Grounding below)

Arrange the power cord and control cords so that they will not create a tripping hazard and where the controls are located to your liking and are conveniently accessible.

Be sure access to plug is not blocked for disconnecting the table from power.


GROUNDING



DANGER

Risk of Electric Shock - Connect this furnishing to a properly grounded outlet only.

This product must be grounded. If it should malfunction or break down, grounding provides a path of least resistance for electrical current to reduce the risk of electric shock. This product is equipped with a cord having an equipment-grounding conductor and a grounding pin. The pin must be plugged into an appropriate outlet that is properly installed and grounded in accordance with all local codes and ordinances. See U.S. sample below.




Improper connection of the equipment-grounding conductor can result in a risk of electric shock. Check with a qualified electrician or service person if you are in doubt as to whether the product is properly grounded. Do not modify the plug provided with the product - if it will not fit the outlet; have a proper outlet installed by a qualified electrician.

ADJUSTING THE BACK REST SECTION




CAUTION

Do not sit on the Backrest section.



TO RAISE THE BACKREST:

Raise the back rest by pressing the lever under the head section and lifting up to the desired height. It is recommended that clients take some of their weight off the back rest while you lift.



TO LOWER THE BACKREST:

Press the lever under the head section and push the back rest down.

DIRECTIONS FOR USE

CONTROL OPERATIONS - ALL TABLES (EXCEPT PT400 & PT400M)




CAUTION


Be sure the self-locking molded rubber feet are flat side down to ensure safety and stability.

The PT tables offer height positioning with electronic ease. Operate the Foot Control and Optional Hand Control as shown below to raise or lower the height of the table.

FOOT CONTROL OPERATION



HAND CONTROL OPERATION (OPTIONAL)




ARTICULATING HEAD REST USE



CAUTION

Do not sit or stand on the articulating head rest section. Do not use Head Rest Section in assisting the client on or off the table.



TO RAISE THE HEAD REST

Grasp the edge of the section and lift up until you have achieved the angle needed. Release and it will lock into place automatically.



TO LOWER THE HEAD REST:

Twist and hold the knob under the head rest and push the head rest down. Release the knob and the head rest will stay in place.

STANDARD SPECIFICATIONS	
Motion 1	Powered 16"-34" (40-86 cm.) Height Range
Motion 2	Manual 0° - 90° Fowler Backrest
Motion 3	Manual 40° to -90° Articulating Head Rest
Foot Control	Height Motion
Table Voltage Options	120V/60Hz or 230V/50Hz
Table Capacity	550 lbs. (250 kg.) (patient and accessories)
Widths	27" (69 cm.), 29" (74 cm.), or 31" (79 cm.)
Length - 3 Section top	39" (99 cm.) Fowler Backrest 23" (28 cm.) Seat Section 10" (25 cm.) Head Rest
Padding	1.75" (4.4 cm.) Firm Response <i>OR</i> 2.5" (6.4 cm.) Comfort Foam
Table Weight	160 lbs. (73 kg)

OPTIONS
Dual Locking Casters
2 Button Hand Control (Height Motion)
Second Foot Control
ACCESSORIES
(1) or (2) Side Rails



Fabrication Enterprises, Inc.
PO Box 1500, White Plains, NY 10602
1-800-431-2830
www.fab-ent.com