

# Upright Bike 4.0U

**SPIRIT**  
powered by *dyaco*



Experience the transformative power of the 4.0U upright bike, meticulously crafted to rejuvenate rehabilitation with unrivaled comfort and precision. Engineered with precision and comfort at its core, plus self-generating power eliminates the hassle of cords.

Conquer rehabilitation goals with exercises tailored to your client's needs:

- Cordless convenience with self-generating power
- 40 levels of resistance
- Built-in cooling fan, water bottle holder, and CSAFE connection
- 11 preset programs and built-in heart rate monitoring
- Fore and aft seat adjustment with high-density foam seat

The upright position fosters proper posture and alignment, maximizing the effectiveness of every session.



Seat can be adjusted vertically and horizontally



Upper handlebars provide added stability during intense workouts



Intuitive, easy-to-read display

## Specifications

Resistance System	Electromagnetic resistance system
Drive Train	Heavy-duty, 8-groove poly V belt with spring-loaded idlers
Power	Self-generating
Workload Range	10 to 750 Watts
Recommended User Height	4'10" - 6'6"
Max User Weight	440 lb
Dimensions	42.5" x 21.6" x 53.5"
Product Weight	133 lb
Compliance	CE Class SA UL (Commercial grade) EN957
Warranty	5 Years parts, 1 Year labor