



FABRICATION  
ENTERPRISES INC

## Model 15-4207 Bariatric Training Stairs

### ***Bariatric Staircase Features:***

- 48" wide access, with Dual Platforms (48" x 48" & 36" x 48").
- 1000 LB capacity, with tubular steel, hardwood and engineered hardwood construction.
- Triple handrail design at 35", 26" and 17" for adults and children.
- Step height 6", tread depth 10", with anti-slip tread and platform surfaces
- Floor Footprint 78" x 104"; Max height at upper platform handrail is 60".
- Minimum ceiling height 102" to 108".
- Ships "knocked down", assembly required
- Ship weight 750 lbs.



**Bariatric Training Stairs**  
*48" Wide - Dual Platforms*

