







STAINLESS STEEL TRAINING STAIRS MODEL TS-200

STANDARD FEATURES

- All Stainless Steel Construction
- 1-1/2" Diameter Handrails
- Steps descend from both sides
- Slip Resistant Material adhered to Tread



TS-200

Fixture May Show Some Available Options

SPECIFICATIONS

Unit is fabricated of 14 gage type 304 stainless steel with exposed exterior surfaces polished to a satin finish. Provides 11" deep Treads with six 4-1/4" Risers on one side and four 6-1/4" Risers on the opposite side.

Model: TS-200
 Top Platform: 31" x 31"
 Overall Dims: 119" L x 34-1/4" W x 59-5/8" H
 Weight Capacity: 800 lbs.

Please visit www.whitehallmfg.com for most current specifications.



Member of
U.S. Green
Building
Council



Recyclable
Product