



Division of
Acorn Engineering Company




A MEMBER OF



STAINLESS STEEL TRAINING STAIRS MODEL TS-400

STANDARD FEATURES

- All Stainless Steel Construction
- 1-1/2" Diameter Handrails
- Steps descend from one side
- Slip Resistant Material adhered to Tread





TS-400

Fixture May Show Some Available Options

SPECIFICATIONS

Unit is fabricated of 14 gage type 304 stainless steel with exposed exterior surfaces polished to a satin finish. Provides 11" deep Treads with four 6-1/4" Risers on one side.

Model: TS-400
 Top Platform: 31" x 31"
 Overall Dims: 64" L x 34-1/4" W x 59-5/8" H
 Weight Capacity: 800 lbs.



Please visit www.whitehallmfg.com
for most current specifications.



Member of
U.S. Green
Building
Council



Recyclable
Product