



## Evolv Large Packages

### Evolv Large User Specifications

Approximate Height Range	5'0"-6'2"
Weight Limit	280 lbs.
Seat Depth Range (seat to back surface)	18"-23"
Seat to Foot Plate Range (from seat pivot)	14"-21"

#### Minimum Support Package

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Removable Flat Back, Velcro® Positioning Belt, Chest Strap, Black Molded Tray with Chest Pad, Manual Hydraulic Actuator with Handle, Flip-Up Knee Pad, Planar Seat, Multi-Adjustable Foot Plates, 5" Front Wheels, 5" Rear Locking Casters.

#### Actuator Handle Location (select one):

- Right (standard)
- Left

#### Back (select one):

- Removable Flat Back 19"H (standard)
- Removable Flat Back 15"H

#### Positioning Belt (select one):

- Velcro® Positioning Belt 26"-55" (standard)
- Velcro® Positioning Belt 38"-67"

#### Chest Strap (select one):

- Chest Strap 34"-54" (standard)
- Chest Strap 30"-46"
- Chest Strap 46"-66"

#### Moderate Support Package

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Swing-Away Front, Adjustable/Removable Actuator Handle, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads for Swing-Away, Removable Flat Back, Velcro® Positioning Belt, Chest Strap, Accessories Mounting Bracket, Lateral Supports, Planar Seat, Multi-Adjustable Foot Plates, Black Molded Tray, 5" Front Wheels, 5" Rear Locking Casters.

#### Foot Straps (select one):

- Secure Foot Straps 15"L (length over top of foot) (standard)
- Secure Foot Straps 21"L (length over top of foot)

#### Hip Supports (select one):

- Hip Supports 12"-17" range (standard)
- Hip Supports 10"-15" range
- Hip Supports 17"-23" range

#### Back (select one):

- Removable Flat Back 19"H (standard)
- Removable Flat Back 15"H

#### Positioning Belt (select one):

- Velcro® Positioning Belt 26"-55" (standard)
- Velcro® Positioning Belt 38"-67"

#### Chest Strap (select one):

- Chest Strap 34"-54" (standard)
- Chest Strap 30"-46"
- Chest Strap 46"-66"

#### Lateral Supports (select one):

- Lateral Supports 11"-19.5"W (standard)
- Lateral Supports 8"-16.5"W

#### Maximum Support Package

The maximum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Front Swivel Casters for Swing-Away Front, Swing-Away Front for Shadow Tray, Adjustable/Removable Actuator Handle, Shadow Tray, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads for Swing-Away, Removable Contoured Back, Positioning Belt, Accessories Mounting Bracket, Chest Vest, High Mount Chest Vest Bracket, Lateral Supports, Head Support, Black Molded Tray, Multi-Adjustable

#### Shadow Tray (select one):

- Shadow Tray 10" (10"H from seat, 4"-15"D from back) (standard)
- Shadow Tray 7" (7.25"H from seat, 4"-15"D from back)
- Deep Shadow Tray 7" (7.25"H from seat, 8"-19"D from back)
- Deep Shadow Tray 10" (10"H from seat, 8"-19"D from back)
- Shadow Tray 12" (12.25"H from seat, 4"-15"D from back)
- Deep Shadow Tray 12" (12.25"H from seat, 8"-19"D from back)

#### Foot Straps (select one):

- Secure Foot Straps 15"L (length over top of foot) (standard)
- Secure Foot Straps 21"L (length over top of foot)

#### Hip Supports (select one):

- Hip Supports 12"-17" range (limited adjustment with shadow tray) (standard)
- Hip Supports 10"-15" range
- Hip Supports 17"-23" range (limited adjustment with shadow tray)

#### Back (select one):

- Removable Contoured Back 19"H (standard)
- Removable Contoured Back 15"H

#### Positioning Belt (select one):

- Velcro® Positioning Belt 26"-55" (standard)
- Velcro® Positioning Belt 38"-67"

#### Chest Vest (select one):

- X-Style Chest Vest 14.5"Lx10.5"W
- X-Style Chest Vest 11"Lx9.5"W
- Y-Style Chest Vest 17.5"Lx11.5"W
- Y-Style Chest Vest 20.5"Lx12.5"W

#### Lateral Supports (select one):

- Lateral Supports 11"-19.5"W (standard)
- Lateral Supports 8"-16.5"W

#### Head Support (select one):

- Head Support-Tall (standard)
- Head Support-Short

#### Minimum Support Package (Mobile)

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Mobile, Clear Angle Adjustable Tray for Mobile, Large Contoured Chest Pad, Removable Flat Back, Velcro® Positioning Belt, Chest Strap, Manual Hydraulic Actuator with Handle, Planar Seat, Flip-Up Knee Pads, Multi-Adjustable Foot Plates, 5" Rear Locking Casters.

#### Actuator Handle Location (select one):

- Right (standard)
- Left

#### Back (select one):

- Removable Flat Back 19"H (standard)
- Removable Flat Back 15"H

#### Positioning Belt (select one):

- Velcro® Positioning Belt 26"-55" (standard)
- Velcro® Positioning Belt 38"-67"

#### Chest Strap (select one):

- Chest Strap 34"-54" (standard)
- Chest Strap 30"-46"
- Chest Strap 46"-66"