



Glider Medium & Large Packages

Glider Medium User Specifications

Approximate Height Range	4'0"-5'6"
Weight Limit	200 lbs.
Seat Depth Range (seat to back surface)	20"-25"
Seat to Foot Plate Range (from seat pivot)	13.5"-16.5"

EasyStand Glider Medium

Moderate Support Package

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Adjustable/Removable Actuator Handle, Secure Foot Straps, Hip Supports, Glide Handle Extensions, Removable Flat Back, Velcro® Positioning Belt, Chest Strap, Full Range Seat, Positioning Bar, 5" Front Wheels, 5" Rear Locking Casters, Tray & Chest Pad, Flip-Up Knee Pads, Adjustable Resistance Cylinders, Adjustable Foot Plates.

Foot Straps (select one):

- Secure Foot Straps 15"L (length over top of foot, two pair) (standard)
- Secure Foot Straps 21"L (length over top of foot, two pair)

Hip Supports (select one):

- Hip Supports 9"-15" range (standard)
- Hip Supports 11"-17" range

Back (select one):

- Removable Flat Back 19"H (standard)
- Removable Flat Back 15"H

Positioning Belt (select one):

- Velcro® Positioning Belt 26"-55" (standard)
- Velcro® Positioning Belt 38"-67"

Chest Strap (select one):

- Chest Strap 34"-54" (standard)
- Chest Strap 30"-46"
- Chest Strap 46"-66"

Glider Large User Specifications

Approximate Height Range	5'0"-6'2"
Weight Limit	280 lbs.
Seat Depth Range (seat to back surface)	24"-29"
Seat to Foot Plate Range (from seat pivot)	15.5"-18.5"

EasyStand Glider Large

Moderate Support Package

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Adjustable/Removable Actuator Handle, Secure Foot Straps, Hip Supports, Removable Flat Back, Velcro® Positioning Belt, Chest Strap, Full Range Seat, Positioning Bar, 5" Front Wheels, 5" Rear Locking Casters, Tray & Chest Pad, Flip-Up Knee Pads, Adjustable Resistance Cylinders, Glide Handles, Adjustable Foot Plates.

Foot Straps (select one):

- Secure Foot Straps 21"L (length over top of foot, two pair) (standard)
- Secure Foot Straps 15"L (length over top of foot, two pair)

Hip Supports (select one):

- Hip Supports 12"-17" range (standard)
- Hip Supports 10"-15" range

Back (select one):

- Removable Flat Back 19"H (standard)
- Removable Flat Back 15"H

Positioning Belt (select one):

- Velcro® Positioning Belt 38"-67" (standard)
- Velcro® Positioning Belt 26"-55"

Chest Strap (select one):

- Chest Strap 34"-54" (standard)
- Chest Strap 30"-46"
- Chest Strap 46"-66"