



MPS Size 1 Packages

MPS Size 1 User Specifications

Approximate Height Range up to 44"
Weight Limit 70 lbs.

Patterns and Colors

Hygienic Cover Pattern (select one):

- Space Race (standard)
- Rainbow Unicorn
- Sea Creatures
- Fancy Pants Elephants
- All Smiles Animals

Frame Color (select one):

- Green (standard)
- Purple
- White
- Orange

Minimum Support Package

The minimum support package includes: (no substitutions)

Mast, Platform Foot Plate, Foot Holders, Foot Straps, Calf/Knee Pad with Strap, Form to Fit Pelvic Pad, Hygienic Covers for (Form to Fit Pelvic Pad, Form to Fit Upper Body Pad, Head Support), Form to Fit Upper Body Pad, Padded Positioning Strap (2), Head Support, Black Molded Swing Away Tray, Swivel Locking Casters, Foot Operated Gas Spring, Mast Receiver with Prone Position Lockout.

Foot Holders (select one):

- Foot Holders 7.75"Lx3.25"W (standard)
- Foot Holders 6"Lx3"W
- Foot Holders 9.75"Lx4"W

Foot Straps (select one):

- Foot Straps 10"L (length over top of foot) (standard)
- Foot Straps 8"L (length over top of foot)

Pelvic Support (select one):

- Form to Fit Pad 5"Hx10"W (standard)
- Form to Fit Pad 5"Hx12"W

Upper Body Support (select one):

- Form to Fit Pad 5"Hx10"W (standard)
- Form to Fit Pad 5"Hx12"W
- Form to Fit Pad 7"Hx10"W
- Form to Fit Pad 7"Hx12"W

Head Support (select one):

- Head Support-Form to Fit 5"Hx10"W (standard)
- Head Support 5"Hx8"W
- Head Support 6"Hx10"W

Moderate Support Package

The moderate support package includes: (no substitutions)

Mast with Leg Abduction, Multi-Adjustable Foot Plates, Secure Foot Straps, Multi-Adjustable Knee Pads, Multi-Adjustable Calf Pads, Planar Pelvic Pad, Hip Supports, Padded Positioning Strap, Hygienic Covers for (Planar Pelvic Pad, Planar Upper Body Pad, Lateral Supports, Head Support), Planar Upper Body Pad, Lateral Supports, Chest Vest, Head Support, Black Molded Swing Away Tray, Forearm Wings, Push Handle, Swivel Locking Casters, Foot Operated Gas Spring, Mast Receiver with Prone Position Lockout.

Foot Plates (select one):

- Multi-Adjustable Foot Plates 7.75"Lx3.25"W (standard)
- Multi-Adjustable Foot Plates 6"Lx3"W
- Multi-Adjustable Foot Plates 9.75"Lx4"W

Foot Straps (select one):

- Secure Foot Straps 10"L (length over top of foot) (standard)
- Secure Foot Straps 8"L (length over top of foot)

Knee Pads (select one):

- Multi-Adjustable Knee Pads 3.25" (inside width) (standard)
- Multi-Adjustable Knee Pads 2.5" (inside width)
- Multi-Adjustable Knee Pads 4.25" (inside width)
- Multi-Adjustable Knee Pads 5" (inside width)

Calf Pads (select one):

- Multi-Adjustable Calf Pads 3.5" (standard)
- Multi-Adjustable Calf Pads 4.5"

Pelvic Support (select one):

- Planar Pad 5"Hx9"W (standard)
- Planar Pad 5"Hx7"W

Upper Body Support (select one):

- Planar Pad 5"Hx9"W (standard)
- Planar Pad 5"Hx7"W
- Planar Pad 7"Hx7"W
- Planar Pad 7"Hx9"W

(Moderate Support Package Continued)

Lateral Supports (select one):

- Lateral Supports (standard)
- Lateral Supports with Elbow Stop
- Elbow Stop with Arm Rest

Chest Vest (select one):

- X-Style Chest Vest 9.5"Lx9"W (standard)
- X-Style Chest Vest 11"Lx9.5"W

Head Support (select one):

- Head Support-Form to Fit 5"Hx10"W (standard)
- Head Support 5"Hx8"W
- Head Support 6"Hx10"W

Maximum Support Package

The maximum support package includes: (no substitutions)

Mast with Leg Abduction, Dual Control, Gas Spring Lift Lockout, Directional Locking Caster, Multi-Adjustable Foot Plates, Secure Foot Straps, Multi-Adjustable Knee Pads, Multi-Adjustable Calf Pads, Planar Pelvic Pad, Hip Supports, Padded Positioning Strap, Hygienic Covers for (Planar Pelvic Pad, Planar Upper Body Pad, Lateral Supports, Head Support), Planar Upper Body Pad, Lateral Supports, Chest Vest, Head Support, Clear Swing Away Tray, Forearm Wings, Hand Grips, Push Handle, Tool Pouch, Swivel Locking Casters, Foot Operated Gas Spring, Mast Receiver with Prone Position Lockout.

Foot Plates (select one):

- Multi-Adjustable Foot Plates 7.75"Lx3.25"W (standard)
- Multi-Adjustable Foot Plates 6"Lx3"W
- Multi-Adjustable Foot Plates 9.75"Lx4"W

Foot Straps (select one):

- Secure Foot Straps 10"L (length over top of foot) (standard)
- Secure Foot Straps 8"L (length over top of foot)

Knee Pads (select one):

- Multi-Adjustable Knee Pads 3.25" (inside width) (standard)
- Multi-Adjustable Knee Pads 2.5" (inside width)
- Multi-Adjustable Knee Pads 4.25" (inside width)
- Multi-Adjustable Knee Pads 5" (inside width)

Calf Pads (select one):

- Multi-Adjustable Calf Pads 3.5" (standard)
- Multi-Adjustable Calf Pads 4.5"

Pelvic Support (select one):

- Planar Pad 5"Hx9"W (standard)
- Planar Pad 5"Hx7"W

Upper Body Support (select one):

- Planar Pad 5"Hx9"W (standard)
- Planar Pad 5"Hx7"W
- Planar Pad 7"Hx7"W
- Planar Pad 7"Hx9"W

Lateral Supports (select one):

- Lateral Supports (standard)
- Lateral Supports with Elbow Stop
- Elbow Stop with Arm Rest

Chest Vest (select one):

- X-Style Chest Vest 9.5"Lx9"W (standard)
- X-Style Chest Vest 11"Lx9.5"W

Head Support (select one):

- Head Support-Form to Fit 5"Hx10"W (standard)
- Head Support 5"Hx8"W
- Head Support 6"Hx10"W