




*shown with TG0542.EC Tapered Extension

PARALLEL BAR, BARIATRIC



IEC60601-1(ed.3)
RM 14971



1000 lb capacity

Motorized
Height & Width Adjustable

<u>Item No.</u>	<u>Length</u>	<u>Platform Width</u>	<u>Bar Range</u>	<u>Between Uprights</u>
TWG1000.BA	6 ft.	54"W	21 1/2" - 30 1/2"	37 1/2"
TWG1001.BA	8 ft.	54"W	21 1/2" - 30 1/2"	37 1/2"
TWG1002.BA	10 ft.	54"W	21 1/2" - 30 1/2"	37 1/2"
TWG1004.BA	12 ft.	54"W	21 1/2" - 30 1/2"	37 1/2"
TWG1006.BA	16 ft.	54"W	21 1/2" - 30 1/2"	37 1/2"

FEATURES:

- 1000 lb Weight Capacity
- UL Classified 60601-1
- Motorized Height & Width Adjustable
- Fingertip Controls, Individual or Simultaneous Operation
- Tapered Platform Extension Accessories Available
- Steel Frame

Bar Range: Horizontal Width Adjustment: inside dimension, 21 1/2" - 30 1/2"
Vertical Height Adjustment: 25 1/2" - 39 3/4" to the top of the bar

Base: 10 Gauge Steel; 1 3/4" high
12" Tapered End and/or 24" Platform Extension, Optional

Handrails: Anodized Aluminum; 1 5/8" Standard O.D.;

Finish: Pearl Grey Standard; Porcelain Optional. Chromium plated Supports.

Electrical: 115V, 6.8A, 782W, hospital grade plug

Matting: Non-skid black mat: width & length to fit.

Warranty: Ten (10) Year Limited Warranty.

*Nominal Dimensions Stated