

# Tubigrip®

## Multi-purpose tubular support bandage

Provides firm support for sprains, strains and swelling



## What size to choose?

| Tubigrip size | INCHES                            |                                   | CENTIMETERS                     |                                 | Sub-bandage pressure (mmHg) |
|---------------|-----------------------------------|-----------------------------------|---------------------------------|---------------------------------|-----------------------------|
|               | Minimum limb circumference (inch) | Maximum limb circumference (inch) | Minimum limb circumference (cm) | Maximum limb circumference (cm) |                             |
| A             | 3.9                               | 5.7                               | 10                              | 14.5                            | 20 - 30                     |
| B             | 5.7                               | 11.4                              | 14.5                            | 29                              | 10 - 20                     |
| C             | 8.7                               | 13.8                              | 22                              | 35                              |                             |
| D             | 10.2                              | 15.7                              | 26                              | 40                              |                             |
| E             | 11.4                              | 18.1                              | 29                              | 46                              |                             |
| F             | 13.8                              | 20.9                              | 35                              | 53                              |                             |
| G             | 15.7                              | 24.4                              | 40                              | 62                              | 5 - 15                      |
| J             | 21.3                              | 32.3                              | 54                              | 82                              |                             |
| K             | 27.6                              | 44.5                              | 70                              | 113                             |                             |
| L             | 41.7                              | 58.7                              | 106                             | 149                             |                             |
| M             | 46.9                              | 65.7                              | 167                             | 167                             |                             |

\* Tubigrip Tape Measuring Guide available upon request

Tubigrip is available under a full range of sizes to suit most applications. The sub-bandage pressures are obtained in double layer.<sup>1</sup>

Reference: 1. Mölnlycke Health Care. Data on file.

The Mölnlycke and Tubigrip names and respective logos are registered trademarks of Mölnlycke Health Care AB. Distributed by Mölnlycke Health Care US, LLC, Peachtree Corners, Georgia 30092. © 2021. Mölnlycke Health Care AB. 1-800-843-8497 MHC-2021-xxxx

Lit# 270651 04.21

