

Chattanooga® clinical traction tables and units

TX® clinical traction units

- easy to operate, digital touch screen interface
- intermittent, static and cyclic traction programs
- user-defined hold, rest, treatment times
- real-time digital treatment tracking
- built-in safety features:
 - patient interrupt switch
 - cervical maximum poundage warning
 - audible signal at end of treatment and when patient interrupt switch is activated
- traction tension parameters 0 - 200 lb
- progressive and regressive settings
- store up to 10 user protocols and up to 14 patient treatment sessions on patient data card
- patient pain mapping

TX® traction unit

00-4759 TX® traction unit 3,515.00

TX® traction stand

00-3600 mobile traction stand 555.00

color options



TX® with mobile traction stand

TX® data card

TX® display can pivot

TX® ADP-400 hi-lo tables

- hi-lo with hand control and factory installed casters
- optional foot control
- 400 lb capacity
- 77"L x 27"W x 21½" - 41½"H
- package contains table, TX® traction unit, flexion stool (TXS-1), deluxe universal belt system, Saunders® cervical traction system (00-7040K) (page 135)
- specify color



traction tables

ADP400H	ADP-400H	table w/ hand switch	4,620.00
ADP400HC	ADP-400HC	table w/ hand switch and front casters	4,780.00
ADP400F	ADP-400F	table w/ foot switch	4,725.00
ADP400FC	ADP-400FC	table w/ foot switch and front casters	4,890.00

traction table packages

00-6874K	ADP-400H	table and traction package	9,150.00
00-6874KC	ADP-400HC	table, front casters and traction package	9,310.00

traction table options

08-6049	retractable front casters	205.00
06-5988	hi-lo foot switch	102.50
04-8082	patient interrupt switch	140.00
06-0000	table paper holder	75.00

TX® TTFT-200 fixed height table

- 32" fixed height; 350 lb capacity
- package contains table, TX® traction unit, flexion stool (TXS-1), deluxe universal belt system, Saunders® cervical traction system (00-7040K) (page 135)
- specify color



traction tables

00-6840 TTFT-200 table 2,100.00

traction table package

00-6878K TTFT-200 table and traction package 6,640.00

TX® TTET-200 hi-lo table

- electric hi-lo, 400 lb capacity
- 77"L x 28"W x 19" - 38"H
- package contains table, TX® traction unit, flexion stool (TXS-1), deluxe universal belt system, Saunders® cervical traction system (00-7040K) (page 135)
- specify color



traction tables

00-6870H	TTET-200H	table w/ hand switch	3,655.00
00-6870HC	TTET-200HC	table w/ hand switch & casters	3,815.00
00-6870F	TTET-200F	table w/ foot switch	3,765.00
00-6870FC	TTET-200FC	table w/ foot switch & casters	3,925.00

traction table package

00-6876K TTET-200H table and traction package 8,190.00