

treatment tables

Moveo™ XP table

combines the benefits of traditional tilt table protocols with active exercise

- hi-lo 20"-37"
- 0°-30° gradual tilt
- 3-70% body weight bearing
- battery back-up
- sliding pelvic section
- fold away foot plates
- exercise band attachment points
- 100"L x 28"W x up to 37"H
- 400 lb capacity
- graystone upholstery only

00-8080 Moveo™ XP table **MAP** 13,025.00



bariatric hi-lo tables

3-section hi-lo treatment table

head section: 15" x 36" raises 45°, lowers 90° and has contoured face/nose opening

opening center section: 22" x 36" does not raise

foot section: 39" x 36" raises to 70°

- 34 oz heavy-duty vinyl with PermaBlock 3 bacterial protection
- 1¾" firm density foam top
- footswitch control
- 76"L x 36"W x 22"- 38"H
- 800 lb capacity
- specify color from chart below


15-1511 3-section table 3,150.00



electric hi-lo steel mat platform table

- electric hi-lo steel frame
- 2" firm density foam top
- 84"L x 48"W x 20"- 30"H
- 900 lb capacity
- specify color from chart below

15-1554 electric hi-lo table 3,500.00




manual hi-lo steel mat platform table

- adjustable hi-lo steel frame with crank hydraulic
- 2" firm density foam top
- 84"L x 48"W x 20"- 30"H
- 900 lb capacity
- specify color from chart below

15-1558 manual hi-lo table 2,300.00



color options for all above



MAP
(minimum advertised
price applies)