


mat tables

fixed height mat tables and raised rim platform tables

upholstered tables


- 2" polyurethane mat is upholstered with vinyl coated nylon
- mat surface is 21"H to accommodate wheelchair transfer
- 500 lb capacity



15-2000	3 x 6 feet	650.00
15-2001	3 x 7 feet	550.00
15-2002	4 x 7 feet	885.00
15-2003	5 x 7 feet	1,100.00
15-2004	6 x 8 feet	1,315.00

raised-rim platforms

- raised-rim to hold mat in place
- optional 2" mat sold separately
- platform surface 19"H
- 500 lb capacity



15-2030	3 x 6 feet	515.00
15-2031	3 x 7 feet	515.00
15-2032	4 x 7 feet	590.00
15-2033	5 x 7 feet	805.00
15-2034	6 x 8 feet	1,020.00

upholstery color

DEFAULT COLOR

royal blue gray white

tan black kelly green

wood finish

DEFAULT FINISH

natural dark walnut traditional cherry

height adjustable hi-lo tables and mat platforms tables

manual tables

- 2" polyurethane mat is upholstered with vinyl coated nylon
- height adjusts from 19"-27"
- manual crank with hydraulic assist
- 350 lb capacity



15-2016	3 x 7 feet	2,700.00
15-2017	4 x 7 feet	2,800.00
15-2018	5 x 7 feet	2,900.00
15-2019	6 x 8 feet	3,000.00

manual platforms

- raised-rim to hold mat in place
- optional 2" mat sold separately
- height adjusts from 19"-27"
- manual crank with hydraulic assist
- 350 lb capacity



15-2046	3 x 7 feet	2,550.00
15-2047	4 x 7 feet	2,600.00
15-2048	5 x 7 feet	2,700.00
15-2049	6 x 8 feet	2,800.00

economy electric tables

- height adjusts from 19"-27"
- 2" upholstered polyurethane mat
- rocker switch activator is electric with hydraulic assist
- 350 lb capacity



15-2026	3 x 7 feet	4,500.00
15-2027	4 x 7 feet	4,600.00
15-2028	5 x 7 feet	4,700.00
15-2029	6 x 8 feet	4,900.00

economy electric platforms


- height adjusts from 19"-27"
- optional 2" mat sold separately
- rocker switch activator is electric with hydraulic assist
- 350 lb capacity



15-2056	3 x 7 feet	4,425.00
15-2057	4 x 7 feet	4,500.00
15-2058	5 x 7 feet	4,575.00
15-2059	6 x 8 feet	4,650.00

deluxe electric tables

- height adjusts from 23"-32"
- 2" upholstered polyurethane mat
- foot pedal activator is fully electric
- 500 lb capacity



15-2082	4 x 7 feet	6,000.00
15-2083	5 x 7 feet	6,200.00
15-2084	6 x 8 feet	6,400.00


deluxe electric platforms

- height adjusts from 23"-32"
- optional 2" mat sold separately
- foot pedal activator is fully electric
- 500 lb capacity



15-2072	4 x 7 feet	5,800.00
15-2073	5 x 7 feet	5,900.00
15-2074	6 x 8 feet	6,150.00

mats for raised rim platforms



- 2" polyurethane mat is upholstered with vinyl-coated nylon for patient comfort
- fits raised-rim platform tables

15-2100	3 x 6 feet	220.00
15-2101	3 x 7 feet	247.50
15-2102	4 x 7 feet	267.50
15-2103	5 x 7 feet	302.50
15-2104	6 x 8 feet	467.50

wall mounted tables

locked position



extended position



- fold and lock table against the wall to save space
- extends only 7" from the wall
- 2" mat is upholstered with heavy-duty, vinyl-coated nylon
- specify mat color: navy blue, green, white, yellow or black

15-2060	3 x 6 feet, 20" H	850.00
15-2062	4 x 7 feet, 20" H	940.00
15-2063	5 x 7 feet, 20" H	1,075.00