


hardwood treatment tables

plain treatment table

- standard treatment table
- 300 lb capacity



	L x W x H	
15-1010	72" x 24" x 30"	402.50
15-1011	78" x 24" x 30"	455.00
15-1012	72" x 30" x 30"	455.00
15-1013	78" x 30" x 30"	510.00

table with shelf

- modified H-brace with shelf
- 300 lb capacity



	L x W x H	
15-1020	72" x 24" x 30"	550.00
15-1021	78" x 24" x 30"	600.00
15-1022	72" x 30" x 30"	600.00
15-1023	78" x 30" x 30"	650.00

enclosed treatment table

- table with enclosed cabinet, sliding doors and drawer
- choose upholstered or raised rim
- 78"L x 30"W x 30"H
- 350 lb capacity



raised-rim table

15-1083	upholstered table	1,100.00
15-1082	raised rim treatment table	1,150.00
15-2123	mat for raised rim table	292.50

hi-ho treatment tables

manual tables



2-section table with drawer



raised-rim table with drawer

- 78"L x 30"W; adjusts 25-33"H using manual crank

plain top with shelf			raised rim top		
	with drawer	2,600.00		with drawer	2,600.00
	without drawer	2,550.00		without drawer	2,550.00
adjustable back with shelf			mat for raised rim top		
	(2-section top)			78" x 30" x 2"	292.50
	with drawer	2,750.00			
	without drawer	2,650.00			

upholstery color

DEFAULT COLOR



wood finish

DEFAULT FINISH



table with H-brace

- H-brace for added stability
- 300 lb capacity



	L x W x H	
15-1000	72" x 24" x 30"	400.00
15-1001	78" x 24" x 30"	450.00
15-1002	72" x 30" x 30"	450.00
15-1003	78" x 30" x 30"	500.00

adjustable back table

- adjustable back raises to 85°
- modified H-brace with shelf
- 300 lb capacity



15-1030	78"L x 30"W x 30"H	860.00
---------	--------------------	--------

space-saving exam table

- exam table with enclosed cabinet, drawer and sliding doors
- space-saving table has fold down 3-position head and foot rests
- folded: 52"L x 24"W x 30"H
- extended: 72"L
- 300 lb capacity



15-1073	enclosed exam table	1,075.00
---------	---------------------	----------

electric tables



plain top table



2-section table with shelves

- 78"L x 30"W; adjusts from 27-39"H; 2" top
- foot control; 110V, 60Hz, 6.5 amp, 720 watt
- 300 lb capacity

plain top			2-section adjustable top		
	plain top	4,675.00		2-section top	7,925.00
	w/shelves	4,975.00		w/shelves	5,275.00