

work tables

circular hi-lo work table



solid butcher block

- foot operated hydraulic pump raises 1¼" butcher-block top from 28"-36"H
- large area underneath table provides enough room for four wheelchairs

dimensions

15-3300	36" diameter, 23" base	990.00
15-3304	42" diameter, 23" base	1,100.00
15-3301	48" diameter, 27" base	1,210.00
15-3302	60" diameter, 27" base	1,540.00

rectangular hi-lo work table




solid butcher block

- foot operated hydraulic pump raises 1¼" butcher-block top from 28"-36"H
- H-frame base provides stability
- worktable is wheelchair accessible
- stabilizer bar prevents rotation


W x L

15-3310	30" x 48"	1,155.00
15-3311	30" x 60"	1,265.00
15-3313	36" x 48"	1,320.00
15-3314	36" x 60"	1,430.00
15-3316	48" x 60"	1,650.00

work and O.T. tables



O.T. model
(apron cut out)
A




work model
B

- heavy duty work tables; O.T. style tables have apron cutouts for wheelchairs
- 30"H is standard

O.T. (A)	work (B)	W x L	
15-3260	15-3290	30" x 30"	275.00
15-3261	15-3291	30" x 48"	325.00
15-3262	15-3292	30" x 60"	350.00
15-3263	15-3293	30" x 72"	425.00
15-3264	15-3294	30" x 96"	550.00
15-3281	15-3271	36" x 48"	350.00
15-3282	15-3272	36" x 60"	425.00
15-3283	15-3273	36" x 72"	500.00
15-3284	15-3274	36" x 96"	625.00

economy hi-lo work tables



butcher block laminate top



oak laminate top

- accommodates wheelchairs
- manual hydraulic tables adjust 28"-35½"H
- built in floor levelers
- black enamel steel constructed frame

15-3250	English oak laminate 48"W x 60"L x 1¼" thick	2,550.00
15-3251	butcher block laminate 48"W x 96"L x 1½" thick	2,950.00

pedestal hi-lo work tables



solid butcher block
rectangular table

- work table powered by electro-mechanical drive
- operating range: 30"-46"H
- 2¼" maple butcher block top
- heavy-duty black and chrome steel base
- foot control is standard; hand control is optional
- 110V, 60Hz, 4.9 amp, 520 watt

15-3255	rectangular table: 48"W x 60"L x 2¼" thick	4,100.00
15-3256	circular table: 60" diameter (not shown)	3,725.00