

## tilt tables

- tilt tables meet the demands of continuous use in busy therapy departments
- platform is 28"W x 78"L with 1" vinyl covered padding
- metal walk off foot plate, 14" x 17½" has safety tread
- three wide safety straps
- special 4" semi-pneumatic casters, with individual step on brakes
- welded steel construction

### specify upholstery



### economy manual



- manual unit designed for home use
- adjusts from horizontal to vertical
- horizontal height is 32" off ground
- black metal frame
- **blue vinyl standard** or specify color from chart

🚚 15-3040 economy 1,800.00

### professional manual



- tilts from full vertical to 12° below horizontal
- horizontal height is 32" off ground
- adjustment retractable handle
- chrome plated frame
- **black vinyl standard** or specify color from chart

🚚 15-3041 manual 2,750.00

### professional electric



- quiet, electric motor assures smooth, vibration-free adjustment
- adjusts from vertical to 8° below horizontal
- horizontal height is 32" off ground
- pneumatic hand-held control is accessible from either side; 7'L cord
- chrome plated frame
- UL hospital grade plug, 10'L power cord
- **black vinyl standard** or specify color from chart

🚚 15-3042 electric 4,175.00

### professional hi-lo electric



- two independent heavy duty electric motors operate height and tilt
- electrical tilt mechanism and indicator 0-90°
- horizontal height adjusts from 25" to 39"
- single pneumatic hand control provides safety and control for both vertical and horizontal positions
- steel frame
- UL hospital grade plug
- **gray only**

🚚 15-3043 electric hi-lo 8,620.00