

## training stairs

### bariatric stairs

- 1000 lb capacity
- hardwood construction with tubular steel
- triple handrails: 35", 26" and 17"H
- 78" x 104" footprint; dual platform, anti-slip treads
- 102" minimum ceiling height
- assembly required; 750 lb

15-4207 bariatric stairs 4,750.00



### convertible stairs



right angle configuration



straight configuration

- can be converted from a straight line to right angle
- anti-slip treads
- handrails 18½", 27", 35"H
- overall height 54"; height to platform 24"; long section 5'10"
- 500 lb capacity

#### platform / short section

15-4204 24x24" / 4'6" 2,300.00  
15-4201 30x30" / 5' 2,350.00  
15-4205 36x36" / 5'6" 2,425.00

### stairs with platform

- straight staircase with platform enclosed on three sides
- anti-slip treads
- overall: 55"L x 30"W x 54"H  
platform: 30"W x 24"D x 24"H
- 4 steps with 6"-10" risers
- includes bus step
- 500 lb capacity



### straight low stairs

- straight staircase for areas with low ceilings
- 2 and 3 step sides with bus step
- anti-slip treads
- overall: 55"L x 30"W x 40"H; platform: 30"W x 24"D x 12"H
- 2 steps with 6"-10" risers
- 500 lb capacity



### straight high stairs

- straight staircase
- anti-slip treads
- overall: 96"L x 36"W x 60"H  
platform: 36"W x 36"D x 24"H
- 8 steps with 6"-10" risers
- includes bus step
- 500 lb capacity



15-4200 training stairs 1,175.00 15-4202 training stairs 1,125.00 15-4203 training stairs 1,125.00

### training stools

- wheelchair training stools for teaching to get in and out of wheelchairs without assistance
- secure together in any sequence
- top size 20" x 10¾"; height of stools 20", 16", 12", 8", 4"
- 250 lb capacity



15-4221 5-stool training set 425.00

### climbing stools



- nested climbing stools with rubber tread
- stools measure 10" x 12" x 2", 12" x 14" x 4", 14" x 18" x 6", and 16" x 20" x 8"
- 250 lb capacity

15-4220 4-stool climbing set 310.00

### curb and ramp set



- curb and ramp training set
- 3 curbs with handles can be nested and stored under ramp platform
- 44" x 84" platform; ¾"-8" incline; curb sizes 40" x 40" x 6", 36" x 36" x 4", and 32" x 32" x 2"; with anti-slip tread
- 300 lb capacity

15-4210 curb/ramp training set 1,250.00

### stackable stools

- set of 3 stools: 12"W x 24"L each
- non-slip tread, natural hardwood construction
- 3 heights: 4", 6" and 8", stack height 18"
- easy carry handle
- 500 lb capacity

15-4222 stack stools 350.00

