




nosey cups

- nose cutout encourages correct head position
- allows user to drink with little or no shoulder flexion, wrist extension or head/neck movement even when wearing a cervical collar
- dishwasher-safe plastic



A	60-1052	8 oz, sandstone	4.50
B	60-1051	8 oz, light blue	4.50
C	60-1050	8 oz, clear	4.50



■ tint color varies			
D	60-1040	4 oz, clear	6.00
E	60-1041	8 oz, clear	4.50
F	60-1042	12 oz, clear	7.00



G	60-1043	2-handled, 8 oz	8.00
----------	---------	-----------------	------

Kennedy J-cup



- lightweight, spillproof drinking cup
- won't spill even when turned upside down
- ergonomic "J" handle
- use with any straw
- dishwasher safe

60-1000	7 oz cup with lid	9.00
60-1001	replacement lid	9.00

no-tip cup holder



- prevents cups from tipping
- fits cups and glasses with 2 1/4" base
- 5" diameter at base

60-1250	no-tip cup holder	7.75
---------	-------------------	------



spillproof cups and lids



- available in a variety of handle configurations
- each comes with two versatile, specialty lids
- cups and mugs have 8 oz capacity
- useable with or without anti-splash lid
- dishwasher safe

H	60-1010	tumbler	8.00
I	60-1011	1 handle	8.00
J	60-1012	2 handle	8.00

vacuum feeding cup

- permits dribble-free drinking without raising head
- mouthpiece large enough to use with straw
- graduated at 4, 6, and 8 oz
- dishwasher and autoclave safe



60-1210	vacuum cup	3.75
60-1210-10	box of 10	35.00



anti-splash lid

- helps direct liquid into mouth
- can be used with a straw

spout lid

- fits comfortably in mouth
- flow of liquid can be regulated by placing a finger over hole
- dishwasher safe

60-1031	anti-splash	10.00
60-1030	spout	10.00
60-1032	vacuum (not shown)	16.00

halo cup

- wide base provides stability and lessens accidental spills
- unique handle eliminates the need for gripping power
- 6 1/2 oz; stackable
- dishwasher safe



60-1062	halo cup	6.00
---------	----------	------