

hook and loop fastener

hook or loop rolls standard non-adhesive webbing back



- standard woven-back tape is sewn or glued onto object
- hook and loop sold separately

hook	loop	specify white (W), black (BLK), tan (T)	
24-7066	24-7067	10 yd	1/2" x 10 yd 11.25
24-7002	24-7003		1" x 10 yd 15.00
24-7006	24-7007		1 1/2" x 10 yd 18.75
24-7022	24-7023		2" x 10 yd 22.50
24-7064	24-7065	25 yd	1/2" x 25 yd 22.50
24-7000	24-7001		1" x 25 yd 30.00
24-7004	24-7005		1 1/2" x 25 yd 37.50
24-7020	24-7021		2" x 25 yd 45.00

elastic (stretchable) loop rolls



- secure an object with this stretchable elastic loop
- loop can be applied to standard and adhesive-backed hook

specify white (W), black (BLK)		
24-7017	1" x 10 yd	30.00
24-7015	1" x 25 yd	75.00
24-7018	2" x 10 yd	55.00
24-7016	2" x 25 yd	135.00


double-sided loop rolls



- soft, durable, non-stretchable strapping
- loop fabric is on both sides of webbing and either side can mate with any hook

24-7255	1" x 10 yd	30.00
24-7256	1" x 25 yd	47.50
24-7257	2" x 10 yd	47.50
24-7258	2" x 25 yd	95.00

hook circles and tabs (adhesive)




- adhesive-backed tabs or circles
- use to attach splint straps or other applications that require small pieces of hook

rectangular		
24-7245	2"x1" (175/pack)	22.50
circular		
24-7246	3/8" dia (850/pack)	42.50
24-7247	1 1/8" dia (600/pack)	60.00
24-7248	1 1/4" dia (450/pack)	80.00

adhesive back hook or loop rolls


- self-adhesive backing is peeled off and tape is pressed onto an object
- hook and loop sold separately



hook	loop	specify white (W), black (BLK), tan (T)	
24-7056	24-7057	10 yd	1/2" x 10 yd 15.00
24-7012	24-7013		1" x 10 yd 20.00
24-7026	24-7027		1 1/2" x 10 yd 25.00
24-7032	24-7033		2" x 10 yd 30.00
24-7054	24-7055	25 yd	1/2" x 25 yd 30.00
24-7010	24-7011		1" x 25 yd 40.00
24-7024	24-7025		1 1/2" x 25 yd 50.00
24-7030	24-7031		2" x 25 yd 60.00

back-to-back hook and loop rolls


- hook and loop are bonded back-to-back on single tape
- create circular closures for positioning, restraining and strapping



specify white (W), black (BLK), tan (T); white is default color

24-7070	1/2" x 10 yd	15.00	24-7074	1/2" x 25 yd	30.00
24-7071	1" x 10 yd	20.00	24-7075	1" x 25 yd	40.00
24-7072	1 1/2" x 10 yd	30.00	24-7076	1 1/2" x 25 yd	60.00
24-7073	2" x 10 yd	40.00	24-7077	2" x 25 yd	80.00

loop strap with adhesive hook



- regular loop with 2" of adhesive hook on one end

pack of 10		
24-7265	1" x 12"	12.50
24-7266	1" x 18"	15.00
24-7267	1" x 24"	17.50
24-7268	2" x 12"	20.00
24-7269	2" x 18"	22.50
24-7270	2" x 24"	25.00

D-ring loop strap with adhesive hook



- D-ring hook and loop strap
- with or without adhesive material

regular hook pack of 10		
24-7100	1" x 12"	17.50
24-7101	1" x 18"	20.00
24-7102	1" x 24"	22.50
24-7103	2" x 12"	25.00
24-7104	2" x 18"	27.50
24-7105	2" x 24"	30.00
adhesive back		
24-7106	1" x 12"	25.00
24-7107	1" x 18"	27.50
24-7108	1" x 24"	30.00
24-7109	2" x 12"	32.50
24-7110	2" x 18"	35.00
24-7111	2" x 24"	37.50

specify white (W), black (BLK), tan (T); white is default color